

**ROMA**  
MEDICAL



Roma Solva

**ROMA**  
MEDICAL

# ROMA SOLVA



Height Adjustable Seat with Flip Up Armrests



Dismantled for Storage & Transportation.



Removable Battery Pack for Off Board Charging



Easy to Use Controls & Basket

## Special Features

- Detachable battery pack for off board charging.
- Quickly dismantles without tools for easy transport and storage.
- Non marking puncture proof tyres.
- Comfortable adjustable height seat.
- Flip up and width adjustable armrests.
- Ample legroom.
- Available in Red & Blue



## Specifications

Model No.*	<b>R100</b>	Battery Size 12v x 2	<b>12Ah</b>
Class (Invalid Carriage)	<b>2</b>	Range**	<b>12 km / 8 miles</b>
Overall Length	<b>1060 mm / 42"</b>	Weight	<b>54 kgs / 119 lbs</b>
Overall Width	<b>540 mm / 21"</b>	Weight of Heaviest Part (from chassis assembly)	<b>17.5 kgs / 38.5 lbs</b>
Overall Height	<b>860 mm / 34"</b>	Weight of Battery Pack	<b>8.5 kgs / 19 lbs</b>
Seat Width	<b>410 mm / 16"</b>	Weight of Seat	<b>11.5 kgs / 26 lbs</b>
Seat Depth	<b>400 mm / 16"</b>	Maximum Speed	<b>6 km/h / 4 mph</b>
Seat Height (to ground)	<b>560mm-610mm / 22"-24"</b>	Wheel / Tyre Size (Front)	<b>180 mm / 7"</b>
Height Folded (no seat)	<b>490 mm / 19"</b>	(Back)	<b>200 mm / 8"</b>
Backrest Height (from Seat)	<b>370 mm / 14.5"</b>	Turning Radius	<b>1040 mm / 41"</b>
Seat Height (to floorpan)	<b>440-490mm / 17"-19"</b>	Battery Charger (off board)	<b>1.8A</b>
Safe Working Slope	<b>6° / 10%</b>	Ground Clearance	<b>65 mm / 2.5"</b>
Maximum User Weight	<b>136 kg / 21 st</b>		

\* Lower limb amputees should refer to user manual before use. \*\* Passenger weight, ground conditions, temperature and battery condition/size can affect distance of travel.



# ROMA MEDICAL

York Road Bridgend Industrial Estate Bridgend CF31 3TB  
Telephone : 01656 674488 Facsimile : 01656 674499  
Email : sales@romamedical.co.uk

**www.romamedical.co.uk**

Dealer Stamp